Recommended Product

MS® ExcePure Classic Series UF Membrane Filter

This product use anti-pollution of polyinylidenefluoride(PVDF), High chemical stability and strong oxidation resistance, Membrane noodles is fingered hole-anisomerous structure, high anti-pollution, not easy to in-depth pollution. Low pressure operation, low energy consumption, long service life.



MS® High Efficiency Active Carbon Filter

This depth filter were developed to satisfy general purpose carbon adsorption applications in pharmaceutical industry. Free of carbon dust, Simplified handling and cleaning, High economic efficiency due to a long service life.






Validation Service Center

Validation Service

MS has established a country wide network of filtration validation laboratory that are staffed with highly trained and experienced scientists and engineers. MS maintains close working relationships with government and industry organizations around the world to ensure that we are up-to-date with the latest GMP requirements and industry best practices.

MS provide comprehensive knowledge and expertise to examine unique product and process characteristics and to select the correct test methods and conditions to validate your process.

Process of Validation Service





Mambrane Solutions


Membrane Solutions For Pharmaceutical Industry

Pharmaceutical industry

Professional Technical to Filtration Solutions
MS Provide Optimization Solutions to Customers



Formulation LVP/SVP Filtration



Recommended Product

Long Service Life PES Filter Cartridge

MS® hydrophilic PES filter cartridge design by asymmetric hydrophilic PES membrane which has excellent throughput and high durability.High particle removal efficency and economical cost make it widely uesed in pharmaceutical and beverage filtration.




High Efficieny Absolute PTFE Filter Cartridge

MS® FlorPure-GS PTFE filter cartridge has broad chemical compatibility, high particle retention efficiency and high flow rate. Suit for organic solvents, aggress acids, chemical liquid and gas filtration, widely uesed in pharmaceutical industry.



Serum Filtration



Recommended Product

Hydrophilic PVDF Filter Cartridge

MS® hydrophilic PVDF filter cartridge has uniform pore size distribution to ensure effective filtration, high flow rate and low pressure drop.Unique low protein adsorption make it widely used in biopharmaceutical industry.



Positively-charged Nylon66 Filter Cartridge


MS® Positively-charged Nylon66 sterilizing grade filter provides enhanced retention of particles smaller than the pore size (virus, endotoxin),Broad chemical compatibility,High filtration area,combined with low pressure drop,ensure filters high flow rates.




Address: JP Morgan International Plaza Suite 650, 14241 Dallas Parkway, Dallas, TX, US Toll free: 1-214-775-2940 Fax: 1-732-412-4040

Email: info@membrane-solutions.com www.membrane-solutions.com

API (Antibiotics) Filtration



Opthalmic Filtration



Recommended Product

MS® AFS-Pure PFA Fliter Cartridge

MS® AFS-Pure PFA filter cartridge using expanded PTFE membranes and PFA core and cage, will meet all kinds of critical filtration requirement in the most aggressive environments. Its advanced fluid dynamics design provide enhanced lateral flow across, through and out of the membrane.

